

For Research Use Only

# RGS18 Polyclonal antibody

Catalog Number: 11866-1-AP

1 Publications



## Basic Information

Catalog Number:

11866-1-AP

Concentration:

240 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2445

GenBank Accession Number:

BC020632

GeneID (NCBI):

64407

UNIPROT ID:

Q9NS28

Full Name:

regulator of G-protein signaling 18

Calculated MW:

235 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

mouse

Positive Controls:

WB : human peripheral blood platelets,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yong Meng	38214344	Int J Mol Med	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human peripheral blood leukocyte sample were subjected to SDS PAGE followed by western blot with 11866-1-AP (RGS18 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.