

For Research Use Only

# RGS18 Polyclonal antibody, PBS Only

Catalog Number: 11866-1-PBS



## Basic Information

Catalog Number:

11866-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2445

GenBank Accession Number:

BC020632

GeneID (NCBI):

64407

UNIPROT ID:

Q9NS28

Full Name:

regulator of G-protein signaling 18

Calculated MW:

235 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human peripheral blood leukocyte sample were subjected to SDS PAGE followed by western blot with 11866-1-AP (RGS18 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11866-1-PBS in a different storage buffer formulation.