

For Research Use Only

MS4A6A Polyclonal antibody

Catalog Number: 11930-1-AP



Basic Information

Catalog Number:

11930-1-AP

Size:

650 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2464

GenBank Accession Number:

BC022854

GeneID (NCBI):

64231

UNIPROT ID:

Q9H2W1

Full Name:

membrane-spanning 4-domains,
subfamily A, member 6A

Calculated MW:

248 aa, 27 kDa

Observed MW:

18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : K-562 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells lysates were subjected to SDS PAGE followed by western blot with 11930-1-AP (MS4A6A antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.