

For Research Use Only

TRIM6 Polyclonal antibody, PBS Only

Catalog Number: 11953-1-PBS

Featured Product



Basic Information

Catalog Number:

11953-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2556

GenBank Accession Number:

BC065575

GeneID (NCBI):

117854

UNIPROT ID:

Q9C030

Full Name:

tripartite motif-containing 6

Calculated MW:

488 aa, 56 kDa

Observed MW:

56 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cell lysates were subjected to SDS PAGE followed by western blot with 11953-1-AP (TRIM6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11953-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer using 11953-1-AP (TRIM6 antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 11953-1-PBS in a different storage buffer formulation.