

For Research Use Only

SPANXC Polyclonal antibody

Catalog Number: 11957-1-AP



Basic Information

Catalog Number:

11957-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2561

GenBank Accession Number:

BC054023

GeneID (NCBI):

64663

UNIPROT ID:

Q9NY87

Full Name:

SPANX family, member C

Calculated MW:

97 aa, 11 kDa

Observed MW:

20 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human testis tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 11957-1-AP (SPANXC antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.