

For Research Use Only

DTX3L Polyclonal antibody

Catalog Number: 11963-1-AP

3 Publications



Basic Information

Catalog Number:

11963-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2568

GenBank Accession Number:

BC060509

GeneID (NCBI):

151636

UNIPROT ID:

Q8TDB6

Full Name:

deltex 3-like (Drosophila)

Calculated MW:

740 aa, 84 kDa

Observed MW:

84 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, goat

Positive Controls:

WB : MCF-7 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xupeng Zang	33810054	Biomolecules	IHC
Pan Hu	38042777	Cancer Cell Int	WB
Xin-Yu Ding	31854042	Cancer Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 11963-1-AP (DTX3L antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.