

For Research Use Only

SNF8 Polyclonal antibody

Catalog Number: 11980-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

11980-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2606

GenBank Accession Number:

BC008976

GeneID (NCBI):

11267

UNIPROT ID:

Q96H20

Full Name:

SNF8, ESCRT-II complex subunit, homolog (S. cerevisiae)

Calculated MW:

258 aa, 29 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Kattareeya Kumthip	29084253	PLoS Pathog	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

