

For Research Use Only

OLFM2 Polyclonal antibody

Catalog Number: 11997-1-AP



Basic Information

Catalog Number:

11997-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2642

GenBank Accession Number:

BC011361

GeneID (NCBI):

93145

UNIPROT ID:

O95897

Full Name:

olfactomedin 2

Calculated MW:

454 aa, 51 kDa

Observed MW:

51 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human brain tissue were subjected to SDS PAGE followed by western blot with 11997-1-AP (OLFM2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.