

For Research Use Only

GGT5 Polyclonal antibody

Catalog Number: 12002-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

12002-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2641

GenBank Accession Number:

BC011362

GeneID (NCBI):

2687

UNIPROT ID:

P36269

Full Name:

gamma-glutamyltransferase 5

Calculated MW:

586 aa, 62 kDa

Observed MW:

80 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HEK-293 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xing-Cheng Zhao	29996846	BMC Med Genomics	WB
Zhuang Luo	38581025	BMC Med Genomics	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 12002-1-AP (GGT5 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.