

For Research Use Only

TM7SF2 Polyclonal antibody, PBS Only

Catalog Number:12033-1-PBS



Basic Information

Catalog Number:

12033-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2637

GenBank Accession Number:

BC012857

GeneID (NCBI):

7108

UNIPROT ID:

O76062

Full Name:

transmembrane 7 superfamily member 2

Calculated MW:

418 aa, 46 kDa

Observed MW:

55-60 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




human lung tissue were subjected to SDS PAGE followed by western blot with 12033-1-AP (TM7SF2 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12033-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 12033-1-AP (TM7SF2 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12033-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 12033-1-AP (TM7SF2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12033-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed PC-3 cells using TM7SF2 antibody (12033-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 12033-1-PBS in a different storage buffer formulation.