

For Research Use Only

APIP Polyclonal antibody

Catalog Number: 12089-2-AP



Basic Information

Catalog Number:

12089-2-AP

Size:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2729

GenBank Accession Number:

BC009077

GeneID (NCBI):

51074

UNIPROT ID:

Q96GX9

Full Name:

APAF1 interacting protein

Calculated MW:

242 aa, 27 kDa

Observed MW:

27 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : K-562 cells, PC-3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12089-2-AP (AIP antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.