

For Research Use Only

RFX5 Polyclonal antibody

Catalog Number: 12137-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

12137-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2783

GenBank Accession Number:

BC017471

GeneID (NCBI):

5993

UNIPROT ID:

P48382

Full Name:

regulatory factor X, 5 (influences HLA class II expression)

Calculated MW:

616 aa, 65 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : HeLa cells, Jurkat cells, K-562 cells, Raji cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Dong-Bo Chen	32883983	Sci Rep	WB
Dong-Bo Chen	31188160	Chin Med J (Engl)	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12137-1-AP (RFX5 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.