

For Research Use Only

ATL1 Polyclonal antibody

Catalog Number: 12149-1-AP **1 Publications**



Basic Information

Catalog Number: 12149-1-AP	GenBank Accession Number: BC010708	Purification Method: Antigen affinity purification
Size: 130 µg/ml	GeneID (NCBI): 51062	
Source: Rabbit	UNIPROT ID: Q8WXF7	
Isotype: IgG	Full Name: atlastin GTPase 1	
Immunogen Catalog Number: AG2797	Calculated MW: 64 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB, ColP

Species Specificity:
human, mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yupeng Zhao	32916628	iScience	WB, ColP

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data