

For Research Use Only

# EZH2 Polyclonal antibody

Catalog Number: 12151-1-AP **2 Publications**



## Basic Information

Catalog Number:

12151-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2799

GenBank Accession Number:

BC010858

GeneID (NCBI):

2146

UNIPROT ID:

Q15910

Full Name:

enhancer of zeste homolog 2  
(Drosophila)

Calculated MW:

751 aa, 86 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Wu Yang	40562234	J Control Release	
Tao Zhang	37633478	Cell Signal	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data