

For Research Use Only

GPS2 Polyclonal antibody

Catalog Number: 12177-1-AP

1 Publications



Basic Information

Catalog Number:

12177-1-AP

Size:

490 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2821

GenBank Accession Number:

BC013652

GeneID (NCBI):

2874

UNIPROT ID:

Q13227

Full Name:

G protein pathway suppressor 2

Calculated MW:

327 aa, 37 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ting Hu	38755610	Respir Res	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

