

For Research Use Only

# NLRP1 Polyclonal antibody, PBS Only

Catalog Number: 12256-1-PBS



## Basic Information

Catalog Number:

12256-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2908

GenBank Accession Number:

BC016965

GeneID (NCBI):

22861

UNIPROT ID:

Q9C000

Full Name:

NLR family, pyrin domain containing

1

Calculated MW:

155 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 12256-1-AP (NLRP1 antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12256-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lymphoma tissue slide using 12256-1-AP (NLRP1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12256-1-PBS in a different storage buffer formulation.