

For Research Use Only

INPP1 Polyclonal antibody

Catalog Number: 12261-1-AP



Basic Information

Catalog Number:

12261-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2910

GenBank Accession Number:

BC015496

GeneID (NCBI):

3628

UNIPROT ID:

P49441

Full Name:

inositol polyphosphate-1-phosphatase

Calculated MW:

399 aa, 44 kDa

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12261-1-AP (INPP1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.