

For Research Use Only

# NARG2 Polyclonal antibody, PBS Only

Catalog Number:12272-1-PBS



## Basic Information

Catalog Number:

12272-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2914

GenBank Accession Number:

BC013684

GeneID (NCBI):

79664

UNIPROT ID:

Q659A1

Full Name:

NMDA receptor regulated 2

Calculated MW:

982 aa, 109 kDa

Observed MW:

140 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IF/ICC, IP, Indirect ELISA

Species Specificity:

human

## Background Information

NARG2 (NMDA Receptor Regulated 2), also known as C16orf35, is a protein-coding gene located on human chromosome 16p13.3. It is widely expressed across various tissues, with highest expression in brain, testis, and heart.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



IP result of anti-NARG2 (IP:12272-1-AP, 4ug; Detection:12272-1-AP 1:600) with HEK-293T cells lysate 2080 ug. This data was developed using the same antibody clone with 12272-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using NARG2 antibody (12272-1-AP) at dilution of 1:2000 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 12272-1-PBS in a different storage buffer formulation.