

For Research Use Only

FNTA Polyclonal antibody

Catalog Number: 12274-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

12274-1-AP

Concentration:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2918

GenBank Accession Number:

BC017029

GeneID (NCBI):

2339

UNIPROT ID:

P49354

Full Name:

farnesyltransferase, CAAX box, alpha

Calculated MW:

44 kDa

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : JAR cells, NIH/3T3 cells, 3T3-L1 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Linke Huang	35300293	Front Pharmacol	WB
Zengpeng Chi	37853939	Cell Biol Int	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



JAR cells were subjected to SDS PAGE followed by western blot with 12274-1-AP (FNTA antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.