

For Research Use Only

RAB7B Polyclonal antibody

Catalog Number: 12283-1-AP

2 Publications



Basic Information

Catalog Number:

12283-1-AP

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2936

GenBank Accession Number:

BC017092

GeneID (NCBI):

338382

UNIPROT ID:

Q96AH8

Full Name:

RAB7B, member RAS oncogene family

Calculated MW:

199 aa, 23 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Cited Species:

human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ryohei Yamada	38604071	Biochem Biophys Res Commun	IF
Chuang Li	37961430	bioRxiv	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

