

For Research Use Only

RAB7B Polyclonal antibody

Catalog Number: 12283-1-AP **1 Publications**



Basic Information

Catalog Number: 12283-1-AP	GenBank Accession Number: BC017092	Purification Method: Antigen affinity purification
Size: 250 µg/ml	GeneID (NCBI): 338382	
Source: Rabbit	UNIPROT ID: Q96AH8	
Isotype: IgG	Full Name: RAB7B, member RAS oncogene family	
Immunogen Catalog Number: AG2936	Calculated MW: 199 aa, 23 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB, IF

Species Specificity:
human

Cited Species:
mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Chuang Li	37961430	bioRxiv	WB,IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data