

For Research Use Only

# SEMA4A Polyclonal antibody

Catalog Number: 12288-2-AP **1 Publications**



## Basic Information

Catalog Number:

12288-2-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2981

GenBank Accession Number:

BC020974

GeneID (NCBI):

64218

UNIPROT ID:

Q9H3S1

Full Name:

sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4A

Calculated MW:

761 aa, 84 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Hyun Ju Lee	29129916	Nat Commun	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data