

For Research Use Only

# GATAD2A Polyclonal antibody

Catalog Number: 12294-1-AP

2 Publications



## Basic Information

Catalog Number:

12294-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2947

GenBank Accession Number:

BC012902

GeneID (NCBI):

54815

UNIPROT ID:

Q86YP4

Full Name:

GATA zinc finger domain containing 2A

Calculated MW:

633 aa, 68 kDa

Observed MW:

68 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, ChIP

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : K-562 cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yao Yao	31324198	J Exp Clin Cancer Res	WB
Xue Chen	38676947	Nucleic Acids Res	ChIP

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 12294-1-AP (GATAD2A antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.