

For Research Use Only

ARSG Polyclonal antibody

Catalog Number: 12328-1-AP



Basic Information

Catalog Number:

12328-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2991

GenBank Accession Number:

BC012375

GeneID (NCBI):

22901

UNIPROT ID:

Q96EG1

Full Name:

arylsulfatase G

Calculated MW:

527 aa, 57 kDa

Observed MW:

62 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, HEK-293 cells, mouse cerebellum tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12328-1-AP (ARSG antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with 12328-1-AP (ARSG antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.