

For Research Use Only

QTRT1 Polyclonal antibody

Catalog Number:12349-1-AP

Featured Product



Basic Information

Catalog Number: 12349-1-AP	GenBank Accession Number: BC015350	Purification Method: Antigen affinity purification
Size: 1100 µg/ml	GeneID (NCBI): 81890	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9BXRO	
Isotype: IgG	Full Name: queuine tRNA-ribosyltransferase 1	
Immunogen Catalog Number: AG3010	Calculated MW: 403 aa, 44 kDa	
	Observed MW: 44 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of QTRT1 antibody (12349-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-QTRT1 transfected HeLa cells.