

For Research Use Only

EBI3/IL-27B Polyclonal antibody, PBS Only

Catalog Number: 12371-1-PBS

Featured Product



Basic Information

Catalog Number:

12371-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3051

GenBank Accession Number:

BC015364

GeneID (NCBI):

10148

UNIPROT ID:

Q14213

Full Name:

Epstein-Barr virus induced 3

Calculated MW:

229 aa, 25 kDa

Observed MW:

25-35 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IP, Indirect ELISA

Species Specificity:

human, pig

Background Information

EBI3, a 34-kDa secreted glycoprotein, is a subunit in 2 distinct heterodimeric cytokines: IL-27 and IL-35. IL-27, composed of p28 (IL-27) and EBI3, can trigger signaling in T cells, B cells, and myeloid cells. IL-35, composed of EBI3 and the ubiquitously expressed p35 subunit of IL-12, is an inhibitory cytokine involved in regulatory T-cell function. EBI3 plays a role in regulating cell-mediated immune responses. (PMID: 8551575; 18033300; 17874423; 21931777)

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data




human placenta tissue was subjected to SDS PAGE followed by western blot with 12371-1-AP (EBI3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12371-1-PBS in a different storage buffer formulation.



IP result of anti-EBI3 (IP:12371-1-AP, 4ug; Detection:12371-1-AP 1:800) with human placenta tissue lysate 600 ug. This data was developed using the same antibody clone with 12371-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human spleen tissue slide using 12371-1-AP (EBI3 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12371-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human spleen tissue slide using 12371-1-AP (EBI3 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12371-1-PBS in a different storage buffer formulation.