

For Research Use Only

RAB34 Polyclonal antibody

Catalog Number: 12378-1-AP

1 Publications



Basic Information

Catalog Number:

12378-1-AP

Size:

290 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3036

GenBank Accession Number:

BC016841

GeneID (NCBI):

83871

UNIPROT ID:

Q9BZG1

Full Name:

RAB34, member RAS oncogene family

Calculated MW:

251 aa, 28 kDa

Observed MW:

29 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Pavlos Nathan J NJ	20926670	J Neurosci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

