

For Research Use Only

IL-17RA Polyclonal antibody

Catalog Number: 12402-1-AP



Basic Information

Catalog Number:

12402-1-AP

Size:

170 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2955

GenBank Accession Number:

BC011624

GeneID (NCBI):

23765

UNIPROT ID:

Q96F46

Full Name:

interleukin 17 receptor A

Calculated MW:

866 aa, 96 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

IL17RA, also named as CD217 and IL17R, is a receptor for IL17A and IL17F. It is a ubiquitously expressed receptor that is essential for IL-17 biologic activity.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

