

For Research Use Only

# DPM1 Polyclonal antibody, PBS Only

Catalog Number: 12403-2-PBS

Featured Product



## Basic Information

Catalog Number:

12403-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1128

GenBank Accession Number:

BC007073

GeneID (NCBI):

8813

UNIPROT ID:

O60762

Full Name:

dolichyl-phosphate  
mannosyltransferase polypeptide 1,  
catalytic subunit

Calculated MW:

260 aa, 30 kDa

Observed MW:

27-32 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

## Background Information

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 12403-2-AP (DPM1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12403-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 12403-2-AP (DPM1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12403-2-PBS in a different storage buffer formulation.