

For Research Use Only

SOCS6 Polyclonal antibody

Catalog Number: 12425-1-AP **2 Publications**



Basic Information

Catalog Number: 12425-1-AP	GenBank Accession Number: BC020082	Purification Method: Antigen affinity purification
Size: 150 µg/ml	GeneID (NCBI): 9306	
Source: Rabbit	UNIPROT ID: O14544	
Isotype: IgG	Full Name: suppressor of cytokine signaling 6	
Immunogen Catalog Number: AG3096	Calculated MW: 60 kDa	
	Observed MW: 60 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Liang Chen	27510094	Oncol Rep	
Guoqing Wang	36623431	Neurotoxicology	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data