

# FLRT1 Polyclonal antibody

Catalog Number: 12440-1-AP

## Basic Information

**Catalog Number:**

12440-1-AP

**Size:**

130 µg/ml

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG3081

**GenBank Accession Number:**

BC018370

**GeneID (NCBI):**

23769

**UNIPROT ID:**

Q9NZU1

**Full Name:**

fibronectin leucine rich  
transmembrane protein 1

**Calculated MW:**

674 aa, 74 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human, mouse, rat

## Background Information

This gene encodes a member of the fibronectin leucine rich transmembrane protein (FLRT) family, functioning in cell adhesion and/or receptor signalling. The encoded protein, expressed in kidney and brain, migrated as a 90-kD protein on SDS-polyacrylamide gels.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

