

For Research Use Only

STK11IP Polyclonal antibody

Catalog Number: 12464-1-AP



Basic Information

Catalog Number:

12464-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3134

GenBank Accession Number:

BC034051

GeneID (NCBI):

114790

UNIPROT ID:

Q8N1F8

Full Name:

serine/threonine kinase 11
interacting protein

Calculated MW:

1099 aa, 121 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung cancer using 12464-1-AP (STK11P antibody) at dilution of 1:50 (under 10x lens).