

For Research Use Only

STK11IP Polyclonal antibody, PBS Only

Catalog Number: 12464-1-PBS



Basic Information

Catalog Number:

12464-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3134

GenBank Accession Number:

BC034051

GeneID (NCBI):

114790

UNIPROT ID:

Q8N1F8

Full Name:

serine/threonine kinase 11
interacting protein

Calculated MW:

1099 aa, 121 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data