

For Research Use Only

RAB40B Polyclonal antibody

Catalog Number: 12471-1-AP **1 Publications**



Basic Information

Catalog Number: 12471-1-AP	GenBank Accession Number: BC018039	Purification Method: Antigen affinity purification
Size: 550 µg/ml	GeneID (NCBI): 10966	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q12829	
Isotype: IgG	Full Name: RAB40B, member RAS oncogene family	
Immunogen Catalog Number: AG3155	Calculated MW: 278 aa, 31 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human gliomas tissue,
Cited Applications: IHC	
Species Specificity: human	
Cited Species: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Martin Zacharias	29976164	BMC Cancer	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human gliomas using 12471-1-AP (RAB40B antibody) at dilution of 1:50 (under 10x lens).