

For Research Use Only

TBL2 Polyclonal antibody

Catalog Number: 12488-1-AP

Featured Product

5 Publications



Basic Information

Catalog Number:

12488-1-AP

Size:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3200

GenBank Accession Number:

BC012938

GeneID (NCBI):

26608

UNIPROT ID:

Q9Y4P3

Full Name:

transducin (beta)-like 2

Calculated MW:

446 aa, 50 kDa

Observed MW:

45-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC 1:20-1:200

IF/ICC 1:10-1:100

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human, mouse

Cited Species:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : COLO 320 cells, HeLa cells, human placenta tissue, Jurkat cells, mouse pancreas tissue, mouse placenta tissue

IP : COLO 320 cells,

IHC : human pancreas cancer tissue, human heart tissue, human pancreas tissue

IF/ICC : HepG2 cells,

Background Information

TBL2 belongs to a member of the beta-transducin protein family. This gene is commonly deleted in Williams-Beuren syndrome, a developmental disorder with multi-system manifestations caused by haploinsufficiency for contiguous genes deleted in chromosome region 7q11.23. TBL2 encodes a protein with four putative WD40 repeats, which may be involved in some intracellular signaling pathway.

Notable Publications

Author	Pubmed ID	Journal	Application
Yoshinori Tsukumo	25393282	PLoS One	WB, IF
Yoshinori Tsukumo	25976671	Biochem Biophys Res Commun	WB, IF
Yoshinori Tsukumo	26239904	J Cell Biochem	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




COLO 320 cells were subjected to SDS PAGE followed by western blot with 12488-1-AP (TBL2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human pancreas cancer using 12488-1-AP (TBL2 antibody) at dilution of 1:50 (under 10x lens).



Immunofluorescent analysis of () fixed HepG2 cells using 12488-1-AP (TBL2 antibody) at dilution of 1:25.



IP result of anti-TBL2 (IP:12488-1-AP, 4ug; Detection:12488-1-AP 1:800) with COLO 320 cells lysate 2400ug.