

For Research Use Only

MAGEA2 Polyclonal antibody

Catalog Number: 12489-1-AP **1 Publications**



Basic Information

Catalog Number: 12489-1-AP	GenBank Accession Number: BC013098	Purification Method: Antigen affinity purification
Size: 160 µg/ml	GeneID (NCBI): 4101	
Source: Rabbit	UNIPROT ID: P43356	
Isotype: IgG	Full Name: melanoma antigen family A, 2	
Immunogen Catalog Number: AG3141	Calculated MW: 314 aa, 35 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Hao Liu	35190529	Cell Death Dis	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data