

For Research Use Only

SGEF Polyclonal antibody, PBS Only

Catalog Number:12493-1-PBS



Basic Information

Catalog Number:

12493-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3198

GenBank Accession Number:

BC016628

GeneID (NCBI):

26084

UNIPROT ID:

Q96DR7

Full Name:

Src homology 3 domain-containing
guanine nucleotide exchange factor

Calculated MW:

97 kDa

Observed MW:

97 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 12493-1-AP (SGEF antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12493-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 12493-1-AP (SGEF antibody) at dilution of 1:200 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12493-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 12493-1-AP (SGEF antibody) at dilution of 1:200 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12493-1-PBS in a different storage buffer formulation.