

For Research Use Only

IL-T3 Polyclonal antibody

Catalog Number: 12524-1-AP

1 Publications



Basic Information

Catalog Number:

12524-1-AP

Size:

170 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3219

GenBank Accession Number:

BC026309

GeneID (NCBI):

11006

UNIPROT ID:

Q8NHJ6

Full Name:

leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 4

Calculated MW:

448 aa, 49 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF, CoIP

Species Specificity:

human

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Heng Zhou	32280997	J Mol Med (Berl)	WB,IF,CoIP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

