

For Research Use Only

PER3 Polyclonal antibody

Catalog Number: 12550-1-AP **5 Publications**



Basic Information

Catalog Number: 12550-1-AP	GenBank Accession Number: BC026102	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 8863	Recommended Dilutions: WB 1:500-1:1000 IHC 1:20-1:200 IF 1:50-1:500
Source: Rabbit	UNIPROT ID: P56645	
Isotype: IgG	Full Name: period homolog 3 (Drosophila)	
Immunogen Catalog Number: AG3400	Calculated MW: 1201 aa, 132 kDa	
	Observed MW: 132 kDa	

Applications

Tested Applications: IF/ICC, IHC, WB, ELISA	Positive Controls: WB : HeLa cells, IHC : human colon cancer tissue, IF : HeLa cells,
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: human, mouse, rat	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

The suprachiasmatic nuclei (SCN) of the anterior hypothalamus function as a circadian clock in the mammalian brain, generating circadian rhythms in physiology and behavior. The period homolog genes *Per3* is an important component of the circadian clock system. *PER3* plays a role in cell proliferation, apoptosis, DNA damage response, cell cycle control, and treatment response of chemotherapy agents in cancers.

Notable Publications

Author	Pubmed ID	Journal	Application
Wenwen Luo	27373683	Int J Biochem Cell Biol	WB
Xiaopeng Song	35765090	Arthritis Res Ther	WB
Xiaomeng Shen	30100061	Biochem Biophys Res Commun	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 12550-1-AP (PER3 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human colon cancer using 12550-1-AP (PER3 antibody) at dilution of 1:50 (under 10x lens).



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using 12550-1-AP (PER3 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).