

For Research Use Only

CLDND1 Polyclonal antibody

Catalog Number: 12567-1-AP



Basic Information

Catalog Number:

12567-1-AP

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3211

GenBank Accession Number:

BC022551

GeneID (NCBI):

56650

UNIPROT ID:

Q9NY35

Full Name:

claudin domain containing 1

Calculated MW:

253 aa, 29 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Claudins are a family of proteins that are the most important components of the tight junctions, where they establish the paracellular barrier that controls the flow of molecules in the intercellular space between the cells of an epithelium. 23 claudins have been identified. They are small (20-27 kilodalton (kDa)) proteins with very similar structure. They have four transmembrane domains, with the N-terminus and the C-terminus in the cytoplasm.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

