

For Research Use Only

THAP1 Polyclonal antibody

Catalog Number: 12584-1-AP

Featured Product

11 Publications



Basic Information

Catalog Number:

12584-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3285

GenBank Accession Number:

BC021721

GeneID (NCBI):

55145

UNIPROT ID:

Q9NVV9

Full Name:

THAP domain containing, apoptosis associated protein 1

Calculated MW:

213 aa, 25 kDa

Observed MW:

25-30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IHC 1:20-1:200

Applications

Tested Applications:

IHC, WB, ELISA

Cited Applications:

WB, IHC, CHIP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: HepG2 cells, human brain tissue, HEK-293 cells, HeLa cells, human heart tissue, A375 cells, SGC-7901 cells, mouse brain tissue, mouse heart tissue

IHC: human colon cancer tissue,

Background Information

THAP1 belongs to the THAP1 family. It is a DNA-binding transcription regulator that regulates endothelial cell proliferation and G1/S cell-cycle progression. THAP1 may also have pro-apoptotic activity by potentiating both serum-withdrawal and TNF-induced apoptosis. Mutations in THAP1 have been associated with dystonia 6. THAP1 encodes a transcription factor with mostly unknown targets. It regulates the expression of DYT1 (TOR1A), another dystonia-causing gene. After characterization of the TOR1A promoter, THAP1 binds to the core promoter of TOR1A. Wild type THAP1 represses the expression of TOR1A, whereas dystonia 6-associated mutant THAP1 results in decreased repression of TOR1A. Catalog# 12584-1-AP is a rabbit polyclonal antibody raised against full-length of human THAP1.

Notable Publications

Author	Pubmed ID	Journal	Application
Maitane Ortiz-Virumbrales	25231164	Acta Neuropathol Commun	WB
Marta Ruiz	26376866	Hum Mol Genet	WB
Fubo Cheng	35015830	Brain	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 12584-1-AP (THAP1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



WB result of THAP1 antibody (12584-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-THAP1 transfected A549 cells.



Immunohistochemical analysis of paraffin-embedded human colon cancer using 12584-1-AP (THAP1 antibody) at dilution of 1:50 (under 40x lens).