

PPP2R1B Polyclonal antibody

Catalog Number: 12621-1-AP

Featured Product

3 Publications

Basic Information

Catalog Number:

12621-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3307

GenBank Accession Number:

BC027596

GeneID (NCBI):

5519

UNIPROT ID:

P30154

Full Name:

protein phosphatase 2 (formerly 2A),
regulatory subunit A, beta isoform

Calculated MW:

667 aa, 74 kDa

Observed MW:

60-74 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IHC 1:50-1:500

IF/ICC 1:50-1:500

Applications

Tested Applications:

WB, IHC, IF/ICC, ELISA

Cited Applications:

WB, IHC, CoIP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

**Note-IHC: suggested antigen retrieval with
TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Positive Controls:

WB : A2780 cells, human kidney tissue, mouse testis
tissue, rat testis tissue, HeLa cells

IHC : human liver tissue,

IF/ICC : Caco-2 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Wen-Ting Hu	36305233	EMBO Rep	WB
Baoyu Chen	32903816	Front Cell Dev Biol	WB
Wei Liu	38429738	Cancer Cell Int	WB,IHC,CoIP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12621-1-AP (PPP2R1B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 12621-1-AP (PPP2R1B antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Ethanol) fixed Caco-2 cells using PPP2R1B antibody (12621-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).