

For Research Use Only

PABPC3 Polyclonal antibody, PBS Only

Catalog Number: 12625-2-PBS

Featured Product



Basic Information

Catalog Number:

12625-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3317

GenBank Accession Number:

BC027617

GeneID (NCBI):

5042

UNIPROT ID:

Q9H361

Full Name:

poly(A) binding protein, cytoplasmic
3

Calculated MW:

631 aa, 70 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12625-2-AP (PABPC3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12625-2-PBS in a different storage buffer formulation.



IP result of anti-PABPC3 (IP:12625-2-AP, 4ug; Detection:12625-2-AP 1:300) with mouse testis tissue lysate 3600ug. This data was developed using the same antibody clone with 12625-2-PBS in a different storage buffer formulation.