

For Research Use Only

DOCK2 Polyclonal antibody

Catalog Number: 12654-1-AP



Basic Information

Catalog Number: 12654-1-AP	GenBank Accession Number: BC016996	Purification Method: Antigen affinity purification
Size: 130 µg/ml	GeneID (NCBI): 1794	
Source: Rabbit	UNIPROT ID: Q92608	
Isotype: IgG	Full Name: dedicator of cytokinesis 2	
Immunogen Catalog Number: AG3246	Calculated MW: 212 kDa, 38 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data