

For Research Use Only

TULP2 Polyclonal antibody

Catalog Number: 12668-1-AP



Basic Information

Catalog Number:

12668-1-AP

Size:

160 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3350

GenBank Accession Number:

BC026070

GeneID (NCBI):

7288

UNIPROT ID:

O00295

Full Name:

tubby like protein 2

Calculated MW:

520 aa, 59 kDa

Observed MW:

55-59 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PC-3 cells, human liver tissue, human lung tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 12668-1-AP (TULP2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.