

For Research Use Only

# SALL2 Polyclonal antibody

Catalog Number: 12679-1-AP **2 Publications**



## Basic Information

Catalog Number:

12679-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3367

GenBank Accession Number:

BC024245

GeneID (NCBI):

6297

UNIPROT ID:

Q9Y467

Full Name:

sal-like 2 (Drosophila)

Calculated MW:

1007 aa, 105 kDa

Observed MW:

105 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HEK-293 cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Cheng Sun	29593314	Nat Commun	WB
Akihiro Fujikawa	28717188	Sci Rep	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 12679-1-AP (SALL2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.