

For Research Use Only

SHKBP1 Polyclonal antibody

Catalog Number: 12703-1-AP



Basic Information

Catalog Number:

12703-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3386

GenBank Accession Number:

BC022855

GeneID (NCBI):

92799

UNIPROT ID:

Q8TBC3

Full Name:

SH3KBP1 binding protein 1

Calculated MW:

707 aa, 76 kDa

Observed MW:

76 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IF 1:10-1:100

Applications

Tested Applications:

IF/ICC, WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : A431 cells,

IF : HepG2 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 12703-1-AP (SHKBP1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of HepG2 cells, using SHKBP1 antibody 12703-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).