

For Research Use Only

RIN3 Polyclonal antibody

Catalog Number: 12709-1-AP **1 Publications**



Basic Information

Catalog Number: 12709-1-AP	GenBank Accession Number: BC025248	Purification Method: Antigen affinity purification
Size: 700 µg/ml	GeneID (NCBI): 79890	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: Q8TB24	
Isotype: IgG	Full Name: Ras and Rab interactor 3	
Immunogen Catalog Number: AG3390	Calculated MW: 985 aa, 108 kDa	
	Observed MW: 108 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells,
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Valentina Latina	37298634	Int J Mol Sci	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 12709-1-AP (RIN3 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.