

For Research Use Only

Synaptotagmin-3 Polyclonal antibody



Catalog Number: 12726-1-AP

Basic Information

Catalog Number:

12726-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3425

GenBank Accession Number:

BC028379

GeneID (NCBI):

84258

UNIPROT ID:

Q9BQG1

Full Name:

synaptotagmin III

Calculated MW:

590 aa, 63 kDa

Observed MW:

63 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse thymus tissue, SH-SY5Y cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 12726-1-AP (Synaptotagmin-3 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 12726-1-AP (Synaptotagmin-3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.