

For Research Use Only

# LIN7B Polyclonal antibody

Catalog Number: 12732-1-AP



## Basic Information

Catalog Number:

12732-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3437

GenBank Accession Number:

BC027618

GeneID (NCBI):

64130

UNIPROT ID:

Q9HAP6

Full Name:

lin-7 homolog B (C. elegans)

Calculated MW:

207 aa, 23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human heart tissue, A549 cells, human brain tissue, K-562 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human heart tissue were subjected to SDS PAGE followed by western blot with 12732-1-AP (LIN7B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.